


## Specificity of LOEWE® Tobamovirus Antisera

ELISA Cat. No.	07184	07188	07042	07047	07096	07176	07175	07065	07022	07167
Pathogen	Tobamo-screen I	Tobamo-screen II	TMV	ToMV	PMMoV	TMGMV	ToBRFV	PaMMV	CGMMV	ZGMMV
TMV	+++	-	+++	+	-	-	+	-	-	-
ToMV	+++	-	-	+++	-	(+)	(+)	-	-	-
PMMoV	+++	-	-	-	+++	+	-	-	-	-
BePMV	++	-	-	-	-	-	-	-	-	-
TMGMV	+	-	-	-	(+)	+++	-	-	-	-
ToMMV	+	-	-	++	-	-	(+)	-	-	-
ToBRFV	+	-	-	(+)	-	-	+++	-	-	-
PaMMV	-	-	-	-	-	(+)	-	+++	-	-
RMV	+	+++	-	-	-	-	-	-	-	-
TVCV	+	+++	+	(+)	-	-	-	-	-	-
YoMV	-	+++	-	-	-	-	-	-	-	-
SFBV	-	+++	-	-	-	-	-	-	-	-
CGMMV	-	-	-	-	-	-	-	-	+++	-
ZGMMV	-	-	-	-	-	-	-	-	-	+++
KGMMV	-	-	-	-	-	-	-	-	-	-
OdRSV	++	-	-	-	-	-	-	-	-	-

## Specificity of LOEWE® Tospovirus Antisera

ELISA Cat. No.	07510	07501	07504	07505	07508
Pathogen	Tosposcreen	TSWV	TCSV	INSV	IYSV
TSWV	+++	+++	(+)	-	-
TCSV	+++	+	+++	-	-
GRSV	+++	+	+	-	-
CSNV	++	-	(+)	-	-
INSV	+++	-	-	+++	-
CACV	(+)	-	-	-	-
WSMoV	(+)	-	-	-	-
IYSV	-	-	-	-	+++



OD*	< 0.1	0.2 - 0.3	0.3 - 1.0	1.0 - 2.5	2.5 - 3.0
Coding	no reaction	slight cross reaction	medium cross reaction	cross reaction	main reaction

\*determined after one hour of substrate incubation